







## Pump Settings

Select:

**Options > My Pump > Personal Profiles > Pump Settings**

**Quick Bolus**  Off  On

**Increment Type**  Units of Insulin  Grams of Carbohydrates

**Increment Amount**  0.5u  1u  2u  5u  
 2g  5g  10g  15g

**Max Bolus**  units

**Basal Limit**  u/hr

## Control-IQ

Select:

**Options > My Pump > Control-IQ**

**Control-IQ**  Off  On

To turn on Control-IQ, you need to enter your weight & Total Daily Insulin. The weight entered must be no less than 24.9kg or 55 lbs and the Total Daily Insulin must be between 10 - 100 units/day.

**Weight**  kg  lbs

kg  lbs

**Total Daily Insulin**  units

## Sleep Activity Schedules (Control-IQ Technology Only)

Select:

**Options > Activity > Sleep Schedules**

Sleep Schedules	On/Off	Selected Days	Start Time	End Time
1				
2				

## Pump Reminder Settings

Select:

**Options > My Pump > Alerts & Reminders > Pump Reminders**

**Low BG**  Off  On

**High BG**  Off  On

**Remind Me Below**  mmol/L

**Remind Me Above**  mmol/L

**Remind Me After**  mins/hrs

**Remind Me After**  mins/hrs

**After Bolus BG** Off  On

Remind Me After

mins/hrs

**Missed Meal Bolus**

Reminder	On/Off	Selected Days	Start Time	End Time
1				
2				
3				
4				

**Pump Alert Settings**

Select:

**Options > My Pump > Alerts & Reminders > Pump Alerts****Low Insulin** units**Auto-Off** Off  OnTime  hours**Device Settings**

Select:

**Options > Device Settings****Display Settings****Screen Time Out** 15 secs  30 secs  60 secs  120 secs**Time and Date** Time and Date set**Sound Volume****Pump Alarms** High  Medium  Low  Vibrate**Pump Alerts** High  Medium  Low  Vibrate**Reminders** High  Medium  Low  Vibrate**Keypad** High  Medium  Low  Vibrate**Bolus** High  Medium  Low  Vibrate**Quick Bolus** High  Medium  Low  Vibrate**Fill Tubing** High  Medium  Low  Vibrate**CGM Alerts** Vibrate  Soft  Normal  HypoRepeat**Security PIN** Off  On PIN

## My CGM

Select:

### Options > My CGM > CGM Alerts

#### High and Low

**High Alert**  Off  On **Alert Me Above**  mmol/L

**Repeat**  Never  15 mins  30 mins  1 hr  
 2 hrs  3 hrs  4 hrs  5 hrs

**Low Alert**  Off  On **Alert Me Below**  mmol/L

**Repeat**  Never  15 mins  30 mins  1 hr  
 2 hrs  3 hrs  4 hrs  5 hrs

#### Rise and Fall

**Rise Alert**  Off  On

**Rate**  0.11 mmol/L/min  0.17 mmol/L/min

**Fall Alert**  Off  On

**Rate**  0.11 mmol/L/min  0.17 mmol/L/min

**Out of Range**  Off  On **Alert Me After**  mins/hours

Select:

### Options > My CGM > Transmitter ID

**Transmitter ID**

## Site Reminder

Select:

### Options > Load > Site Reminder

**Site Reminder**  Off  On

**Remind Me In**  1 day  2 days  3 days

**Remind Me At**



For more information and support on the t:slim X2 insulin pump, please contact the AMSL Diabetes Technical Support Team on **1300 851 056**.

[amsldiabetes.com.au](https://amsldiabetes.com.au)    



Always read the label and use only as directed. Systems like the t:slim X2 insulin pump with Control-IQ technology are not substitutes for the active management of diabetes. Control-IQ technology should not be used by anyone under the age of six years old. It should also not be used in patients who require less than 10 units of insulin per day or who weigh less than 25 kgs. The Control-IQ feature relies on current Dexcom G6 CGM (sold separately) readings to function and will not be able to predict glucose levels and suspend insulin delivery if a user's CGM is not functioning properly or their pump is unable to receive the CGM signal. Users should always pay attention to their symptoms, actively monitor and manage glucose levels and treat according to their healthcare professionals' recommendations. Fast acting insulin is required for the t:slim X2 insulin pump, consult a healthcare professional to see which product is right for you. The Bluetooth® word mark and logos are registered trade-marks owned by the Bluetooth SIG, Inc. and any use of such marks by Tandem Diabetes Care is under license. All other trademarks and copyrights are the property of their respective owners. AMSL is a subsidiary of Dexcom. ARTG 304681. PR-100-467 January 2023